

| Table 6A | | | |
|-----------------|------------------|--------------------------------|---|
| REF | GTIN Code | REF_DESCRIPTION | Extended description |
| BGC-05 | 8056518500312 | BIO-GEN Cortical Bone Granules | Cortical Granules - 6 btl / 0.5g 0.25-1mm |
| BGC-05s | 8056518500169 | BIO-GEN Cortical Bone Granules | Cortical Granules - 1 btl / 0.5g 0.25-1mm |
| BGM-05 | 8056518500336 | BIO-GEN Mix Granules | Cancellous Cortical Granules - 6 btl / 0.5g 0.25-1mm |
| BGM-05n | 8056518500343 | BIO-GEN Mix Granules | Cancellous Cortical Granules - 1 btl / 0.5g 0.25-1mm |
| BGM-05s | 8056518500176 | BIO-GEN Mix Granules | Cancellous Cortical Granules - 1 btl / 0.5g 0.25-1mm |
| BGM-10 | 8056518500350 | BIO-GEN Mix Granules | Cancellous Cortical Granules - 6 btl / 0.25g 0.25-1mm |
| BGM-100 | 8056518503832 | BIO-GEN Mix Granules | Cancellous Cortical Granules - 6 btl / 1g 0.25-1mm |
| BGM-100s | 8056518500183 | BIO-GEN Mix Granules | Cancellous Cortical Granules - 1 btl / 1g 0.25-1mm |
| BGM-10s | 8056518500206 | BIO-GEN Mix Granules | Cancellous Cortical Granules - 1 btl / 0.25g 0.25-1mm |
| BGM-20 | 8056518500367 | BIO-GEN Mix Granules | Cancellous Cortical Granules - 1 btl / 2g 0.25-1mm |
| BGS-05 | 8056518500459 | BIO-GEN Cancellous Granules | Cancellous Granules - 6 btl / 0.5g 0.25-1mm |
| BGS-05n | 8056518500466 | BIO-GEN Cancellous Granules | Cancellous Granules - 1 btl / 0.5g 0.25-1mm |
| BGS-05s | 8056518500558 | BIO-GEN Cancellous Granules | Cancellous Granules - 1 btl / 0.5g 0.25-1mm |
| BGS-09 | 8056518500473 | BIO-GEN Cancellous Granules | Cancellous Granules - 6 btl / 0.5g 1-2mm |
| BGS-09s | 8056518500596 | BIO-GEN Cancellous Granules | Cancellous Granules - 1 btl / 0.5g 1-2mm |
| BGS-10 | 8056518500480 | BIO-GEN Cancellous Granules | Cancellous Granules - 1 btl / 0.5g 1-2mm |
| BGS-11 | 8056518500497 | BIO-GEN Cancellous Granules | Cancellous Granules - 6 btl / 1g 1-2mm |
| BGS-11s | 8056518500671 | BIO-GEN Cancellous Granules | Cancellous Granules - 1 btl / 1g 1-2mm |
| BGS-20 | 8056518500510 | BIO-GEN Cancellous Granules | Cancellous Granules - 1 btl / 2g 0.25-1mm |
| BGS-22 | 8056518500534 | BIO-GEN Cancellous Granules | Cancellous Granules - 1 btl / 2g 1-2mm. |
| BGS-23 | 8056518500541 | BIO-GEN Cancellous Granules | Cancellous Granules - 6 btl / 1g 2-3mm |
| BGS-23s | 8056518500725 | BIO-GEN Cancellous Granules | Cancellous Granules - 1 btl / 1g 2-3mm |



BIOTECK S.p.A.

Headquarters:

Via E. Fermi, 49 - 36057 Arcugnano (VI) Italy
 Ph: (+39) 0444 289366 - Fax: (+39) 0444 285272
 e-mail: info@bioteck.com - e-mail PEC: certificata@pec.bioteck.com - SDI: SUBM70N
 VAT.n: IT02702750247 - CF: 06857400011 - REA: VI268440 - Sh.Cap. €120.000,00 i.v.

Production Facility:

Via G. Agnelli, 3 - 10020 Riva presso Chieri (TO) Italy
www.bioteck.com

8. Illustrative pictures of falsified devices


Feature



1) Description on bottle label. The falsified device has the wording "RESORBABLE NATURAL BONE" (red arrow). In contrast, the correct wording is shown in the "Original Device" column (red arrows).

Falsified Device



Original Device



| Feature | Falsified Device | Original Device |
|--|---|---|
| <p>2) Reference to NB CE on bottle label. The falsified device has the indication: 0474 (red arrow). The original device has the indication for NB 0477. For batches produced until 2019 NB is 0373.</p> |  |  |

| Feature | | |
|--|---|--|
| <p>3) Cap and bottle closure ring. The falsified device has a light blue cap with grooves along the edges (red arrows). The original device, may have two variants, both with dark blue cap, without groove.</p> | <p><i>Falsified device: front view</i></p>  | <p><i>Original device: front view (2 variants on the market)</i></p>  |
| | <p><i>Falsified device: view from above</i></p>  | <p><i>Original device: view from above (2 variants on the market)</i></p>  |